



noah 4

- ride the wave

Noah Mobile technology - the future of Noah

Noah Mobile technology is a service included with Noah (requires Noah 4.9 or later).

With Noah Mobile technology, mobile apps from HIMSA and HIMSA member companies can easily interact with Noah, just like Windows based modules do today. This means you will soon be able to perform Noah-related tasks on mobile devices such as Android tablets and iPads!



For example, with Noah Mobile technology, you can use an audiometer app on your iPad to:


- Establish a secure connection to Noah, either within your office or at a remote patient site
- Search for and open a patient's data record
- Perform a hearing care test on the chosen patient
- Store the new patient data in the Noah database when you are done

And since each Noah Mobile app communicates directly with Noah System on your PC, there's no need to install Noah on your mobile device.

What kinds of apps can use Noah Mobile technology?

HIMSA member companies are free to decide the best types of Noah Mobile apps to provide and which mobile devices to support.

Potentially these apps can include everything from fitting and measurement apps to journal and questionnaire apps, ordering apps, telemedicine apps and more.



Are there Noah Mobile apps available?

Yes! All HIMSA-certified Noah Mobile apps are listed on the HIMSA website www.himsa.com.

HIMSA member companies are hard at work developing new apps using Noah Mobile technology.

Contact your HIMSA member company for more information on the apps they will provide and the expected release dates.

Will my patient data be secure?

HIMSA uses a secure data transfer process when transmitting data between Noah and a Noah Mobile app. Among other things, this means that:

- your data will never be stored on the Internet.
- you have full control over which apps can access your patient data and how much patient data they can access.



HIMSA