



What's new in Noah System 4?

Over the last few years HIMSA has added several new features to Noah System, making it quicker and easier to work with your patients and patient data.

Here are some of the great reasons to upgrade to the latest Noah System version.

Ready to Support Tele-Health

HIMSA has worked closely with our member companies to create a flexible and secure solution that can quickly evolve to meet changing needs and technologies. The result is the Noah Alert Dashboard.

The screenshot shows the Noah 4 Alert Dashboard interface. The window title is "Noah 4". The dashboard has a header "Alert Dashboard" and two dropdown menus: "Sort by: Most recent" and "From: All alerts". There are two buttons: "Alert Details" (highlighted) and "Session List".

aaa	aaa	Hearing Performance	5/24/2019
aaa	aaa	Hearing Performance	5/24/2019
aaa	aaa	Hearing Performance	5/24/2019
James	La Branche	Hearing Performance	5/13/2019
James	La Branche	Hearing Performance	5/13/2019
Dee	Davies	Hearing Performance	4/2/2019
Peter	Collins	Hearing Performance	3/13/2019
Barbara	DOOLING	Usage Status	2/20/2019
Laura	Palmer	Request for Assistance	2/20/2019
Laura	Palmer	Request for Assistance	2/20/2019
Mary	MCHUGH	Hearing Performance	2/18/2019

The details panel on the right shows:

- Patient : Peter Collins
- Noah User : Unassigned [Assign](#)
- Priority :
- Category : Hearing Performance
- Received : 3/13/2019 2:50:33 AM
- Description : Problems hearing the TV.

At the bottom of the details panel, there are icons for "Journal Module" and a "Remove Alert" button. A "Close" button is located at the bottom right of the dashboard window.

With the Alert Dashboard in Noah System 4.12, you will be able to see all important patient requests from HI manufacturers directly in Noah. You will then be able to sort through the alerts, prioritize them, and solve the most critical problems right away.

Although you can already see the Alert Dashboard in Noah 4.12, you will need to wait a while before you can use it. It will take time for manufacturers to implement and test the communication between their portals and Noah System.



Improved Networking Performance

In Noah 4.10, HIMSA updated the network technology for more stable and robust performance. This means more reliable communication with remote offices, especially those with less stable networks.

Enter the information in the fields below. This information is required for your data security records and will be stored in Noah and the Noah Activity Log.

Intended Recipient:

Purpose of Export:

Save ✓ Cancel ✕

Patient Data Security

In today's world, patient data security is essential. That's why in Noah System 4.9 and 4.10, we added a variety of new tools to help make data security as complete and intuitive as possible. These include:

- Active Directory Integration (Windows) login to better control access to Noah System data
- Encryption option for your Noah database
- The ability to encrypt, password-protect and document Noah exports
- Better control of your Noah activity logs

Fast Data Viewing

HIMSA has added a Fast Data Viewing feature, providing you with immediate access to patient data without opening your Noah modules! By clicking on the eye next to a fitting action for example, you can review instrument settings without having to wait for the fitting module to open! Or click the eye icon to quickly browse the patient's journal entry.

Easily Merge Duplicate Patient Records

Noah System provides a new merge patient feature, allowing you to easily locate and merge duplicate patient records in your Noah database.

Patient Categories

Noah System also includes a comprehensive category feature, allowing you to create and assign user-defined patient categories.

A patient's details screen will clearly show the categories to which the patient belongs. And searching, displaying and exporting patients based on categories is easy.

Noah Mobile Technology

Noah Mobile technology lays a secure ground for communication with Noah compatible apps on mobile devices and in the cloud.

You can already purchase audiometer solutions that take advantage of Noah Mobile technology. And in the future, HIMSA members will use Noah Mobile for questionnaire apps, tele-health applications and more.

Noah 4.12 is a free upgrade for all registered Noah 4 versions. Visit www.himsa.com.

